

**YZ11d Ophthalmoscope**

The portable product has excellent optical performance by improved optical design.

**Features:**

- Improved light path and excellent optical performance.
- Imported halogen bulb provides high and adjustable brightness and good color rendering.
- Smart, portable and convenient for mobile inspection.
- Equipped with special charger, quick and convenient recharging, against over recharging.
- Anti-recharging and anti-discharging function protects the battery.

**Specifications:**

Illumination Style	Large Spot, Small Spot, Slit, Net, Red Free
Diopter Compensation	0 ±1 ±2 ±3 ±4 ±5 ±6 ±8 ±10 ±12 ±16 ±20 -25 -35D
Illumination Source	3.5V/2.8W, Halogen Bulb
Input Voltage	AC220V±22V/50Hz±1Hz、AC110V±11V/60Hz±1Hz
Power(Main Body of Ophthalmoscope)	Li+ Battery Pack 3.7V, Rechargeable

---

**6 6 VISION TECH CO.,LTD.**

Add: 34 Daru Lane, Gusu District, CN-215005 Suzhou, Jiangsu / P.R.China, 98 China

Web Site: [www.66vision.com](http://www.66vision.com)

E-mail: [trade@66vision.com](mailto:trade@66vision.com)