

YZ20P Operation Microscope



The product is smart, compact, portable and applicable in various operations.

Detailed profile

Characteristics:

- Smart and compact body allow easy mobility.
- Three-step magnifications.
- Spring balance system of arm ensures the microscope move freely.
- All the lens are multi-coated, mildewproof and anti-reflective.
- Objectives with different focus and various accessories are optional. It can be applied in ENT, Surgery and Gynaecology.

Specifications:

Eyepiece Magnifications	12.5×
Adjustable Diopter	±5D

6 6 VISION TECH CO.,LTD.

Add: 34 Daru Lane, Gusu District, CN-215005 Suzhou, Jiangsu / P.R.China, 98 China

Web Site: www.66vision.com

E-mail: trade@66vision.com

Adjustable Range for Pupil Distance	50mm~70mm
Objective Lens	f=200
Working Distance	190mm
Magnifications	6X 10X 16X
Illumination Source	12V/100W, cold reflection halogen lamp for medical use
Illumination Type	6° cold light source coaxial illumination
Coaxial Illumination	≥100001x
Reach Radius of Arm	860mm
Adjustable Vertical Range	850mm~1100mm
Voltage	AC220V±10%/ 50Hz±1Hz、AC110V±10%/60Hz±1Hz
Power	120VA
Fuse	AC250V T1.25A 、AC125V T2.5A
Electrical Safety Standard	Executive Standard:IEC601—1,Class I, Type B
Packing Volume	0.283m ³ , 2cartons
Total Gross Weight	61kg

6 6 VISION TECH CO.,LTD.

Add: 34 Daru Lane, Gusu District, CN-215005 Suzhou, Jiangsu / P.R.China, 98 China

 Web Site: www.66vision.com

 E-mail: trade@66vision.com