

## YZ24 Streak Retinoscope



The product can quickly and precisely measure the astigmatism axis. The light can be converged, radiated and paralleled.

### Detailed profile

#### Characteristics:

- The well positioned filament can be rotated for 360° and move upward and downward. Brightness of the streak can be adjusted.
- Quickly and precisely measure the astigmatism axis.
- The light can be converged, radiated and paralleled
- Automatic turning off effectively protects the bulb.
- Two-steps adjustable brightness for bulb: S(Bright)、W(Dark).

#### Specifications:

Working Distance	1m
Streak	3mm-20mm
Streak Rotation	360°
Illumination Source	3V/2.1W, Halogen Bulb
Input Voltage	AC220V $\pm$ 22V/50Hz $\pm$ 1Hz、 AC110V $\pm$ 11V/60Hz $\pm$ 1Hz

---

### 6 6 VISION TECH CO.,LTD.

Add: 34 Daru Lane, Gusu District, CN-215005 Suzhou, Jiangsu / P.R.China, 98 China

Web Site: [www.66vision.com](http://www.66vision.com)

E-mail: [trade@66vision.com](mailto:trade@66vision.com)