

YZ3 Binocular Handheld Slit Lamp Microscope



Detailed profile

Characteristics:

- Converging microscope provides sharp observation.
- Excellent optical performance, total 800g weight, easy to operate and practical enough.

6 6 VISION TECH CO.,LTD.

Add: 34 Daru Lane, Gusu District, CN-215005 Suzhou, Jiangsu / P.R.China, 98 China

Web Site: www.66vision.com

E-mail: trade@66vision.com

Specifications:

Microscope Type	Converging microscope
Objective Magnifications	Two steps magnifications
Magnifications	10X、16X
Diopter Adjustment	±8D
Pupil Distance	50mm~75mm
Working Distance	60mm
Slit Width	0~12mm continuously adjustable
Diameter of Light Spot	0.2mm、1mm、5mm、12mm
Filter	Color Temperature Compensation, Redfree, Cobalt Blue
Slit Angle	±30°
Light Source	LED
Maximal Illumination	≥30000Lx
Power	7.4V 2200mA Li Battery, Rechargeable
Working Time	2.5 hours for common use
Packing Volume	216mm X 100mm X 225mm
Net Weight	900g

6 6 VISION TECH CO.,LTD.

Add: 34 Daru Lane, Gusu District, CN-215005 Suzhou, Jiangsu / P.R.China, 98 China

 Web Site: www.66vision.com

 E-mail: trade@66vision.com