

YZ5J Slit Lamp



This is a new product to fill the optical requirements, which is easy to operate and medium-priced.

Detailed profile

Characteristics:

- Two-step magnifications: 10×、16×.
- All the Optical lenses are moistureproof, mildewproof and anti-reflected treated.
- Excellent optical performances.
- Convenient and compact.

6 6 VISION TECH CO.,LTD.

Add: 34 Daru Lane, Gusu District, CN-215005 Suzhou, Jiangsu / P.R.China, 98 China

Web Site: www.66vision.com

E-mail: trade@66vision.com

Specifications:

Microscope Type	Converging stereoscope microscope
Objective Magnification	Two Steps Magnifications
Eye Piece	10x
Magnifications (Field of View)	10x (ϕ 18mm) 、 16x (ϕ 14.5mm)
Pupil Distance	55mm~82mm
Diopter Adjustment	$\pm 6D$
Slit Width	0mm~9mm adjustable (slit is round when the slit width is 9mm)
Slit Height	1mm~8mm adjustable
Aperture of Diameter	ϕ 9mm、 ϕ 8mm、 ϕ 5mm、 ϕ 3mm、 ϕ 2mm、 ϕ 1mm、 ϕ 0.2mm
Slit Angle	0°~180° adjustable
Slit Inclination	5° , 10° , 15° , 20° , four steps adjustable
Diameter of Light Spot	ϕ 9mm、 ϕ 8mm、 ϕ 5mm、 ϕ 3mm、 ϕ 2mm、 ϕ 1mm、 ϕ 0.2mm
Filter	Heat Absorption, Grey, Redfree, Cobalt Blue
Illumination Bulb	12V/30W Halogen Bulb
Input Voltage	AC220V \pm 22 V ; AC110V \pm 11 V
Electrical Safety Standard	Conform to Standard: IEC601-1, Class I, Type B
Packing Volume	720mm \times 495mm \times 480mm
Total Weight	23Kg
Optional Accessories	R type, T type Applanation Tonometer, Measuring Eye Piece, Trihedron Lens, Gonioscope, Funduscope, Aspheric Lens

6 6 VISION TECH CO.,LTD.

Add: 34 Daru Lane, Gusu District, CN-215005 Suzhou, Jiangsu / P.R.China, 98 China

 Web Site: www.66vision.com

 E-mail: trade@66vision.com