

## YZ5X Slit Lamp Microscope



### Specifications:

Microscope Type	Converging stereoscope
Objective Magnification	Two steps drum magnifications
Eye Piece	10X
Magnifications (Field of View)	10x(Φ 22.5mm) 16x (Φ 14mm)
Pupil Distance	55 mm~82 mm
Diopter Adjustment	± 6D
Slit Width	0mm~14mm adjustable(slit is round when the slit width is 14mm)
Slit Height	1mm~14mm adjustable

Slit Angle	0°~180° adjustable
Diameter of Light Spot	Φ 14mm、Φ 9mm、Φ 5.5mm、Φ 0.3mm
Slit Inclination	None
Filter	Heat Absorption, Redfree, Cobalt Blue
Illumination Bulb	6V20W Halogen Bulb
Input Voltage	AC220V±22 V; AC110V±11 V
Electrical Safety Standard	Conform to Standard IEC601-1,Class I,Type B
Total Weight	23Kg
Packing Volume	720mm´495mm´480mm
Optional Accessories	Applanation Tonometer、Trihedron Lens、Gonioscope、Funduscope、Aspheric Lens