



ACP990L Auto Charts Projector



Chart Types



Accessories:



Features:

1. By the way of double locating, it would correct errors automatically
2. Offering clearer and sharper pictures
3. Increase the picture density; there are 5 icons in one line
4. The charts are bigger: 33*22.5(CM) (in 5M refraction)
5. Projected by wireless remote controller
6. Able to project the charts one by one
7. Standard accessories, the metal board is covered by bias material

Specifications:

Projection distance	2-6M
Projection magnification	30×(in 5M refraction)
Projection size	33*22.5(CM),Φ275MM(in 5M refraction)
Numbers of charts	33
Chart change-over	3 frame/sec
Numbers of masks	1 open,3 horizontal line,8 vertical line,16 single isolation,1 red and green
Tilt range	10°horizontally as well as vertically
Bulb	high-luminance and long-life LED light source
Image shift interval	0.5 second
Automatic shut-off mechanism	5 minutes
Input voltage	AC 110V/60HZ or 220V/50HZ
Power consumption	65W
Dimensions	19.6(W)*29.2(L)*23.4(H)(CM)
Weight	6.5KGS