

# COU500 Ophthalmic Chair and Stand



## Features:

1. This ophthalmic unit with luxurious and graceful design
2. The chair can be driven up and down stably, rotated left and right
3. With Less noise, big loading ability and longer lifespan

## Specifications:

Table size	84*40.5*4(CM)
Table rotation range	90° (stop automatically)
Table movable range	about 33(CM)
Swing arm oscillation angle	up and down 30°
Chair rotation range	60°
Chair up and down range	15(CM)
Chair loading ability	180KGS
Chair backrest lay down pitch	130°
Input voltage	AC220V/AC110V
Lamp	15W, AC220V/AC110V
Fuse	5A, AC150V
Dimensions	140*100*199(CM)
Unit weight	175KGS