

COU980 Ophthalmic Chair and Stand



Features:

1. Luxurious design and graceful structure
2. Y shape table can rotate up and down
3. The table can locate several instruments
4. The control panel can drive the table and the chair easily

Specifications:

Unit size	110*113*220(CM)
Table size	44*55(CM)
Table rotation range	360°
Table up and down range	10(CM)
Swing arm oscillation angle	30°
Chair up and down range	55-66(CM)
Chair loading ability	150KGS
Chair backrest slope	180°
Input Voltage	AC220V/AC110V
Lamp	16W, AC220V/110V
Fuse	5A, AC250V
Unit weight	180KGS