

OPHTHALMOSCOPES

EYE, EAR, NOSE & THROAT

3.5 V Standard Ophthalmoscope

Uses halogen light for true tissue color and consistent, long-lasting illumination. Versatile and reliable at an economical price.

- Halogen light for true tissue color and consistent, long-lasting illumination
- 6 apertures for generalist use: micro, small, and large spot sizes, fixation target, slit aperture, and red-free filter
- 28 focusing lenses with a range of -25 to +40 diopters
- Rubber brow rest prevents scratching of eyeglasses



11710
3.5 V Standard Ophthalmoscope
(head only)

11710 Standard Ophthalmoscope (head only)

3.5 V Standard Ophthalmoscope: Accessories

03000-U 3.5 V Halogen Lamp for 11710 Standard Ophthalmoscope

APERTURES FOR: 3.5 V Standard Ophthalmoscope



Micro-Spot Aperture

allows quick entry into small, undilated pupils.



Fixation Aperture

has graduated crosshairs for measuring eccentric fixation or locating lesions.



Small Aperture

provides excellent view of fundus through an undilated pupil.



Slit

used to determine levels of lesions and tumors.



Large Aperture

standard aperture for a dilated pupil and general examination of the eye.



Red-Free Filter

excludes red rays from examination field for easy identification of veins, arteries and nerve fibers.

Lenses	Filters	Optical System
28 Lenses -25 to +40 Diopters	Red-free	Standard Optics