

## RETINOSCOPE

EYE, EAR, NOSE & THROAT



**18245**  
Elite™ Retinoscope  
Shown with 719  
Lithium Ion Handle

### Elite™ Retinoscope

Makes refractions easier than ever by producing one of the brightest, sharpest, and cleanest streaks available.

- Provides 2.8X the brightness of our standard retinoscopes
- Superior streak quality enables easier and quicker refractions—even in patients with small pupils
- Crossed linear polarizing filter eliminates bothersome glare from trial lenses, helping you perform a more precise exam
- External focusing sleeve improves comfort and maneuverability
- Sealed optics keep the lenses cleaner for more effective exams
- Fits all Welch Allyn 3.5 V power sources
- For those preferring a spot retinoscope, the streak retinoscope converts to a spot retinoscope by simply replacing the streak lamp with a spot lamp

**18245** Elite Streak Retinoscope

**18300** Elite Spot Retinoscope

### Elite Retinoscope: Accessories

**08200-U** Halogen HPX® Streak Replacement Lamp

**08300-U** Halogen HPX Spot Replacement Lamp

**18250** Fixation Cards for dynamic retinoscopy (set of four)



**18250**  
Fixation Cards  
Set of four fixation cards quickly and easily attach to magnetic logo on instrument head for dynamic retinoscopy.