

3.5 V Larynx / Nasopharynx Illuminators

Halogen HPX® illumination and fiber optic light transmission for easy exams.

- Fiber optic light transmission produces cool distal light
- Wide-angle rotatable mirror with pilot aperture
- Illuminators withstand sterilization with steam, water, gas or solutions

27000 3.5 V Complete Larynx Illuminator

27050 3.5 V Complete Nasopharynx Illuminator

3.5 V Nasal Illuminators: Accessories

03100-U 3.5 V Halogen Lamp

27010 3.5 V Larynx Mirror Section Only

27070 3.5 V Fiber Optic Illuminator Only

03100-LED LED Upgrade Kit



27000
Larynx Illuminator