

Fiber Optic Sigmoidoscope Set

Fiber Optic Sigmoidoscopes offer the finest in design and craftsmanship to provide optimal clinical performance.

- Halogen light for true tissue color and consistent, long-lasting illumination
- Fiber optics provide a cool distal ring of light with no reflections, no obstructions
- Made from autoclavable stainless steel for convenience and durability
- Hinged window with a neoprene seal avoids fogging during insufflation
- Graduated tube with removable tip for fiber cleaning
- Compatible with available adapters to fit most Welch Allyn power handles



35303
Fiber Optic Sigmoidoscope Set

35303 Fiber Optic Sigmoidoscope Set

35303 Fiber Optic Sigmoidoscope Set includes the following:

- 32410** 15 mm x 15 cm Fiber Optic Pediatric/Stricture Speculum only with Obturator
- 32820** 19 mm x 25 cm Standard Sigmoidoscope Speculum only with Obturator
- 37023** 23 mm x 7 cm Fiber Optic Anoscope Speculum only with Obturator
- 73211** Light Handle with 3 ft Cord
- 73305** Transformer with 5 ft Cord
- 30200** Insufflation Bulb, Complete (Latex-Free)
- 07800-U** 6 V Halogen Lamp (for use with 73211.6V Rectal Light Handle) (pack of 2)
- 30130** 30 cm Suction Tube (standard)
- 05391-U** Fiber Optic Sigmoidoscope Set (35303) Carrying Case



30200
Insufflation Bulb