

ILLUMINATORS

EYE, EAR, NOSE & THROAT



41100
3.5 V Finnoff Ocular
Transilluminator



43300
3.5 V Curved All-Purpose
Transilluminator

3.5 V Fiber Optic Transilluminators

General fiber optic examination light for checking pupil reflex. Transillumination of frontal and maxillary sinuses.

- Fiber optic light transmission for cool distal light
- Fully sterilizable with steam, water, gas or solutions
- Optional cobalt blue filter

41100 3.5 V Finnoff Ocular Transilluminator

41101 3.5 V Finnoff Ocular Transilluminator with Cobalt Blue Filter

43300 3.5 V Curved All-Purpose Transilluminator

03100-U 3.5 V Halogen Lamp

03100-LED LED Upgrade Kit

3.5 V Halogen Fiber Optic Transilluminators: Accessories

41102 Cobalt Blue Filter only