








5079-135
Adult Stethoscope

Professional Adult Stethoscopes

With its minimal weight and custom-fitting headset and ear tips, the Professional is one of the most comfortable and versatile general-use stethoscopes manufactured today. Its acoustics allow for accurate auscultation of cardiac, pulmonary, and vascular sounds. And the chestpieces in the adult, pediatric, and anesthesiology versions of the Professional enable you to select a stethoscope that best meets your specific needs.

- General-use stethoscope
- Stainless steel double-head (bell and flat diaphragm) chestpiece with nonchill rim
- Rotatable binaurals and interchangeable comfort sealing ear tips
- Latex-free

	5079-135 Adult Stethoscope, 28" (71 cm), Black
	5079-137 Adult Stethoscope, 28" (71 cm), Navy
	5079-289 Adult Stethoscope, 28" (71 cm), Blue
	5079-139 Adult Stethoscope, 28" (71 cm), Burgundy
	5079-285 Adult Stethoscope, 28" (71 cm), Forest Green

Professional Adult Stethoscope: Accessories

5079-123	Adult Diaphragm Disc Only
5079-124	Pediatric Diaphragm Disc Only
5079-180	Adult Bell Nonchill Ring, Black
5079-181	Adult Bell Nonchill Ring, Gray
5079-182	Pediatric Bell Nonchill Ring, Black
5079-183	Pediatric Bell Nonchill Ring, Gray
5079-126	Adult Diaphragm Nonchill Ring, Black
5079-184	Adult Diaphragm Nonchill Ring, Gray
5079-127	Pediatric Diaphragm Nonchill Ring, Black
5079-185	Pediatric Diaphragm Nonchill Ring, Gray
5079-195	Binaural/Spring Assembly and Tubing, 28", Black
5079-196	Binaural/Spring Assembly and Tubing, 28", Green
5079-197	Binaural/Spring Assembly and Tubing, 28", Navy
5079-199	Binaural/Spring Assembly and Tubing, 28", Burgundy

5079-210	Double Head Chestpiece Assembly (Adult), Gray
5079-211	Double Head Chestpiece Assembly (Adult), Black
5079-214	Double Head Chestpiece Assembly (Pediatric), Black
5079-215	Double Head Chestpiece Assembly (Pediatric), Gray
5079-320	Double Head Chestpiece Assembly (Infant), Black
5079-172	Hard Ear Tips (Medium), Pair
5079-173	Hard Ear Tips (Large), Pair
5079-313	Comfort Sealing Ear Tips (Small), Gray
5079-366	Comfort Sealing Ear Tips (Large), Black
5079-170	Comfort Sealing Ear Tips (Large), Gray
5079-367	Comfort Sealing Ear Tips (Large), Navy
5079-368	Comfort Sealing Ear Tips (Large), Burgundy
5079-369	Comfort Sealing Ear Tips (Large), Forest Green