

5079-321
Harvey™ DLX,
Triple Head,
Black










Harvey™ DLX™ Stethoscopes

Both the double- and triple-head Harvey DLX stethoscopes are specifically designed to accentuate important sounds and enable you to hear cardiac, pulmonary, and vascular sounds without relying on difficult variable-pressure techniques. The Harvey DLX is equipped with dual-bore tubing that creates two uninterrupted sound channels and works with the optimally weighted brass chestpieces to transmit a full range of physical sounds. Providing clear, crisp sound, the Harvey DLX helps enhance your auscultation capabilities.

- New longer 28" length full-range, cardiology/pulmonary stethoscope
- Available with double-head (brass bell and flat diaphragm) or rotatable triple-head (brass bell and corrugated diaphragm) chestpieces
- Pediatric chestpieces available
- Chrome-plated brass chestpieces and binaurals
- Rotatable binaurals and interchangeable comfort sealing ear tips
- Triple-leaf binaural spring encased in polyurethane allows size adjustment while preventing breakage
- Latex-free
- Dual lumen



5079-325
Harvey™ DLX, Double
Head, Black

	5079-321	Harvey DLX, Triple Head, Black (28")
	5079-322	Harvey DLX, Triple Head, Burgundy (28")
	5079-323	Harvey DLX, Triple Head, Navy (28")
	5079-324	Harvey DLX, Triple Head, Forest Green (28")
	5079-325	Harvey DLX, Double Head, Black (28")
	5079-326	Harvey DLX, Double Head, Burgundy (28")
	5079-327	Harvey DLX, Double Head, Navy (28")
	5079-328	Harvey DLX, Double Head, Forest Green (28")
	5079-325P	Harvey DLX, Double Head, Black, Pediatric (28")

FLAT DIAPHRAGM

Unique, tightly-contained diaphragm is optimized for higher-frequency, critically important heart sounds like murmurs, clicks, and ejection sounds.



CORRUGATED DIAPHRAGM

The corrugated diaphragm is especially good for detecting lower-frequency sounds, such as rumbles, gallops, and ejection sounds.



TRIPLE-HEAD MODELS

Get complete flexibility and efficiency by placing a bell, flat diaphragm, and corrugated diaphragm at your fingertips.



Harvey™ DLX™ Stethoscopes: Accessories

- 5079-336** Comfort Sealing Ear Tips (Large), Black
- 5079-337** Comfort Sealing Ear Tips (Large), Navy
- 5079-338** Comfort Sealing Ear Tips (Large), Burgundy
- 5079-339** Comfort Sealing Ear Tips (Large), Forest Green
- 5079-113** Hard Ear Tips, Extra Large
- 5079-314** Comfort Sealing Ear Tips (Small), Gray
- 5079-233** Comfort Sealing Ear Tips (Large), Gray
- 5079-108** Name Tag
- 5079-370** Diaphragm Assortment, Harvey Holographic, Adult, Four Colors
- 5079-272** Tubing for 28" Model, Burgundy
- 5079-120** Tubing for 28" Model, Black
- 5079-305** Tubing for 28" Model, Forest Green
- 5079-273** Tubing for 28" Model, Navy
- 5079-274** Binaural/Spring Assembly, Burgundy
- 5079-121** Binaural/Spring Assembly, Black
- 5079-307** Binaural/Spring Assembly, Forest Green
- 5079-275** Binaural/Spring Assembly, Navy
- 5079-104** Flat Diaphragm Only
- 5079-05** Corrugated Diaphragm Disc Only
- 5079-319** Diaphragm Assembly
- 5079-318** Bell Chestpiece Assembly
- 5079-316** Pediatric Chestpiece Assembly
- 5079-315** Harvey DLX Accessory Kit (Includes: Pediatric Chestpiece, Flat Diaphragm, Corrugated Diaphragm, One Pair of Soft Ear Tips, One Pair of Hard Ear Tips)
- 5079-124** Pediatric Diaphragm Disc Only



5079-370
Diaphragm Assortment



5079-120
Black Tubing