

Fiber Optic Laryngoscope Handles

Bright, white HPX® illumination offers a clear color-corrected view.

- Patented HPX technology blends halogen and xenon gases under pressure to produce a bright, white light for true tissue color and consistent, long-lasting illumination
- Up to 40% lighter than our previous design, for increased responsiveness and maneuverability
- Knurled finish ensures secure grip
- Single HPX lamp cartridge design used across all sizes of 2.5 V handles

Lightweight Fiber Optic HPX Handles

60813 2.5 V Medium, uses two C-size batteries

60814 2.5 V Penlight, uses two AA-size batteries

60815 2.5 V Stubby, uses two AA-size batteries

2.5 V Rechargeable Fiber Optic HPX Handle

60713 2.5 V Handle only

3.5 V Rechargeable Fiber Optic HPX Handle

60835 3.5 V Handle only



60813
2.5 V Medium
Handle

60814
2.5 V Penlight
Handle

60815
2.5 V Stubby
Handle



60713
2.5 V Handle only

60835
3.5 V Handle only

| Lightweight Fiber Optic HPX Handles | |
|---|--|
| Part Number | Description |
| 60813 | 2.5 V Medium, uses two C-size batteries |
| 60814 | 2.5 V Penlight, uses two AA-size batteries |
| 60815 | 2.5 V Stubby, uses two AA-size batteries |
| 608125-501 | 2.5 V Halogen HPX Gold Lamp Cartridge Assembly |
| 608107 | Switch Plunger/Lamp Shield |
| 0600-U | 2.5 V Halogen HPX Replacement Lamp |
| 2.5 V Rechargeable Fiber Optic HPX Handle | |
| 60713 | 2.5 V Handle only |
| 71140 | Charger only, 110-130 V, 50/60 Hz |
| 72000 | 2.5 V NiCad Rechargeable Battery |
| 608125-501 | 2.5 V Halogen HPX Gold Lamp Cartridge Assembly |
| 06000-U | 2.5 V Halogen HPX Replacement Lamp |
| 3.5 V Rechargeable Fiber Optic HPX Handle | |
| 60835 | 3.5 V Handle only |
| 71140 | Charger only, 110-130 V, 50/60 Hz |
| 72200 | 3.5 V NiCad Rechargeable Battery |
| 608351-501 | 3.5 V Halogen HPX Gold Lamp Cartridge Assembly |
| 608107 | Switch Plunger/Lamp Shield |
| 00300-U | 3.5 V Halogen HPX Replacement Lamp |

Fiber Optic Laryngoscope Handles: Accessories

| | |
|-------------------|---|
| 608107 | Switch Plunger/Lamp Shield (for use with Lightweight, 2.5 V and 3.5 V Rechargeable Fiber Optic Laryngoscope Handles) |
| 608125-501 | 2.5 V Halogen HPX Silver Lamp Cartridge Assembly (for use with Lightweight and 2.5 V Rechargeable Fiber Optic Laryngoscope Handles) |
| 06000-U | 2.5 V Halogen HPX Replacement Lamp (for use with Lightweight and 2.5 V Rechargeable Fiber Optic Laryngoscope Handles) |
| 72000 | 2.5 V NiCad Rechargeable Battery |
| 72200 | 3.5 V NiCad Rechargeable Battery |
| 00300-U | 3.5 V Halogen HPX Replacement Lamp |
| 71140 | Charger only, 110-130 V, 50/60 Hz |
| 608351-501 | 3.5 V Halogen HPX Gold Lamp Cartridge Assembly |