



MacIntosh Fiber Optic Laryngoscope Blades



E-MacIntosh Fiber Optic Laryngoscope Blades



Miller Halogen Fiber Optic Laryngoscope Blades

## Fiber Optic Laryngoscope Blades

Removable, repolishable light pipes make these an affordable choice in fiber optic laryngoscope blades.

- Single-piece stainless steel construction ensures blade integrity
- Damaged light pathways can be repaired at a fraction of the cost of replacement blades
- Free repolishing service extends product life

**69061** #1 MacIntosh Fiber Optic Laryngoscope Blade

**69062** #2 MacIntosh Fiber Optic Laryngoscope Blade

**69063** #3 MacIntosh Fiber Optic Laryngoscope Blade

**69064** #4 MacIntosh Fiber Optic Laryngoscope Blade

**69211** #1 E-MacIntosh Fiber Optic Laryngoscope Blade

**69212** #2 E-MacIntosh Fiber Optic Laryngoscope Blade

**69213** #3 E-MacIntosh Fiber Optic Laryngoscope Blade

**69214** #4 E-MacIntosh Fiber Optic Laryngoscope Blade

**68065** #00 Miller Halogen Fiber Optic Laryngoscope Blade

**68060** #0 Miller Halogen Fiber Optic Laryngoscope Blade

**68061** #1 Miller Halogen Fiber Optic Laryngoscope Blade

**68062** #2 Miller Halogen Fiber Optic Laryngoscope Blade

**68063** #3 Miller Halogen Fiber Optic Laryngoscope Blade

**68064** #4 Miller Halogen Fiber Optic Laryngoscope Blade

MacIntosh Standard (Lamp) Blades			
Part Number	Blade Size	Overall Length	Inside Length (Base to Tip)
<b>69061</b>	<b>#1</b>	89 mm	63 mm
<b>69062</b>	<b>#2</b>	108 mm	82 mm
<b>69063</b>	<b>#3</b>	126 mm	101 mm
<b>69064</b>	<b>#4</b>	160 mm	135 mm
English MacIntosh Fiber Optic Blades*			
<b>69211</b>	<b>#1</b>	95 mm	70 mm
<b>69212</b>	<b>#2</b>	116 mm	90 mm
<b>69213</b>	<b>#3</b>	135 mm	110 mm
<b>69214</b>	<b>#4</b>	156 mm	130 mm
Miller Halogen Fiber Optic Blades*			
<b>68065</b>	<b>#00</b>	8 mm	36 mm
<b>68060</b>	<b>#0</b>	8 mm	53 mm
<b>68061</b>	<b>#1</b>	8 mm	79 mm
<b>68062</b>	<b>#2</b>	10 mm	132 mm
<b>68063</b>	<b>#3</b>	11 mm	172 mm
<b>68064</b>	<b>#4</b>	13 mm	182 mm

\*Dimensions are approximate